	USER MANUAL	12/12/24	FR
		Page 1 on 27	

**USER MANUAL FOR EXAMINATION
LAMP**

Examination Lamp Led

LAMPE N35 LED

LAMPE LED 9

VERRE & QUARTZ *Technologies*

Phone : +33 1 64 66 20 33– Fax : +33 1 64 77 58 97

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Introduction

Dear customer please read this manual carefully before use the product, to protect yourself and protect the other user from any harm.

This device is a medical device Class I compliant with REGULATION (UE) 2017/745 on medical device (Annex VIII) and its successive amendments and integrations.

The manufacturer declares that this product is in compliance with Annex I (General Safety and Performance Requirement) of regulation (UE) 2017/745 and its subsequent amendments and integration, and certifies this compliance by affixing the CE brand.

The product is classified in risk group 1 according to standard IEC 62471 (photobiological risk).

This user manual applies to following models:

- **LED 9**
- **N35 LED**

Customer service is available for any clarification regarding the product its use identification of spare parts and any question you may have about the device and its operation if you wish to order spare parts and for assistance and assistance and warranty issues.

VERRE & QUARTZ Technologies

Phone. : +33 1 64 66 20 33– Fax : +33 1 64 77 58 97

Email: info@verre-et-quartz.fr – Web: www.verre-et-quartz.com

Espace Vinci, 24 Avenue Graham Bell, 77600 Bussy St Georges France

The content of this notice may be modified by Verre et Quartz Technologies, without prior notice or further obligation, in order to make changes and improvements. Any reproduction or translation in whole or part of this notice is prohibited without the written permission of Verre et Quartz Technologies.

Verre et Quartz Technologies reserves the right to modify, remove or change the data contained in this documentation at any time and for any reason without prior, because Verre et Quartz Technologies is constantly seeking new solutions that lead to the evolution of the product.

Verre et Quartz Technologies reserves the right and duty to provide in terms or shape equipment technology, and performance.

Regarding translations into languages other than French, the French version of this manual shall prevail.

TABLE OF CONTENTS

1	GENERAL INFORMATION	3
1.1	Qualification of personnel	4
1.2	Packing, Transport, Storage, Installation Site Characteristics	4
1.3	Graphical symbols Used on the products and Packaging	5
1.4	EU Declaration of Conformity	6
1.5	Warranty Certificate	7
2	Importance of Personal Safety	8
2.1	Intended Used	9
2.2	Safety Conditions (Side Effect)	9
2.3	Ambient Conditions	10
2.4	Contrôles avant chaque utilisation	10
3	Product installation	10
3.1	Table Installation	11
3.2	Wall Installation	11
3.3	Rail Clamp Installation	11
3.4	Stand Version Installation	11
3.5	Hand piece Installation (For LED 9 Model Only)	12
3.6	First start up	12
3.7	Verification After Installation and Control Operation Before use	
4	Product Description and Operation	12
4.1	Description and Operation N35 LED	13
4.2	Description and Operation LED 9	13
5	Cleaning and Disinfection	14
5.1	Cleaning the Product	14
5.2	Disinfection the product	14
5.3	Sterilization of Handpieces (For LED 9 Only)	15
6	Adjustments	16
6.1	Annual Check by the Operator	16
6.2	Repairs	16
6.3	Clutch Adjustments	16
6.4	Fault Finding	18
6.5	Routine Maintenance	18
6.6	Spare Part List	18
7	Technical Data	18
7.1	Technical data N35	18
7.2	Technical data LED 9	19
8	Circuit Diagram	21
8.1	Circuit Diagram N35	21
8.2	Circuit Diagram LED 9	22
9	Déclaration CEM	23

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 4 on 27	

1 General information

The EM (Electro-Medical) device to which this manual refers is a LAMP for diagnosis or observation.

For the purpose of this description this ME device will be abbreviated in this manual as "Product"

This manual is an integral part of the product, in accordance with REGULATION (UE) 2017/745 and its successive amendments and integration. Read and keep this use and maintenance manual close to the product.

Verre et Quartz Technologies disclaims any responsibility for any damage caused to persons or propriety due to the installation, maintenance and use of the product by unqualified operators.

Qualified operators are those who have attended a course on the installation and use of the product organized by Verre et Quartz Technologies does not authorize third parties to carry out extraordinary maintenance. If a problem is identified, contact Verre et Quartz Technologies.

The installation operations of the product are entirely at the expense and care of the end customer; no cost or liabilities regarding the installation and/or implantation of the product can be attributed to Verre et Quartz Technologies.

The masonry work required to prepare the wall for installing the product and the electrical work needed for the electrical supply to the product must be performed solidly and safety according to industry standards by duly qualified personnel.

AS a non-exhaustive example please note that the following professionals are considered to have the necessary qualifications.

- ⇒ Civil Engineer, (Surveyor Construction company) duly registered with the corresponding professional board (for masonry work)
- ⇒ Electrical Engineer, Qualified Electrician Technician licensed to practice as an electrician (for electrical work)

The product is an EM electro -medical device therefore falling under standard IEC 62353 Consequently any operation carried out on the product must be performed in accordance with CEI 62353, where applicable.

1.1 Qualification of personnel

This section describes the requirements and qualifications that individuals involved in various phase of the product's life and use must possess.

Installation	Installer and/or qualified technician
Use	Professional medical personnel
Routine maintenance	Technician with technical professional skills
Additional maintenance	Verre et Quartz Technologies or authorized dealer
Assistance	Verre et Quartz Technologies or authorized dealer
Cleaning	Professional medical and paramedical personnel
Démolition	Comply with current regulations regarding waste disposal. This product must not be disposed of in regular garbage bins. To avoid harming the environment and health due to the release of pollutants into the air, separate the various internal components such as iron, aluminum plastic, and electrical materials and entrust them to designated centers that will handle the recycling.

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Box containing the product along with the user manual

The transport is carried out by all road haulage companies provided that it meets the following characteristics :



Temperature (°C) : -15 / +60 ; Humidity : 10 / 95 % ; atmospheric pressure (hPa) : 500 / 1060.

The Packaged product must be storage (warehouse) in a drive place that complies with the following characteristics:

Temperature (°C) : +10 / +40 ; Humidity : 30 / 75 % ; atmospheric pressure (hPa) : 700 / 1060.

1.3 **Graphic symbols used on the product and the packaging**

Description on symbols and nameplates on the product and in this manual.

	CE Marking attesting the product's compliance with REGULATION EU 2017/745 and subsequent amendments and integration		IMQ Brand
	Manufacturing date (year /month)		Medical device
	Manufacturer's adress		Model number
	RECYCLING ! The product must be recycled separately		Serial number
	Standby		Swiss Agent
	Functional earth		CLASSE II Equipment
	On		Off

 Connection point for line conductor  Top side  Fragile packaging  Do not stack packaging  Humidity to be respect (max limit top right and min limit bottom left)  General warning signal  Instruction must Be Read	 Connection point for conductor  Packaging weight  Protect from rain  Temperature limit (max limit top right and limit min bottom left)  Pressure to be respected (max limit top right and min limit bottom left)  General compulsory signal
--	---

1.4 Conformity Declaration

Established pursuant to Article 19 and Annex IV of REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 concerning medical devices and amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC

Manufacturer : **Verre et Quartz Technologies P. LONGONI S.r.l.**

Registered office address: Via Monterosa, 18/20/22 - 20831 SEREGNO (MB) - ITALY

Unique Registration Number (SRN): **IT-MF-000009224**

This EU Declaration of Conformity is issued under the sole responsibility of the manufacturer.

BASE BUD-ID: **++B880LUMINAIREPM**

Name and commercial designation of the product : **ALFA-FIX, ALFA-FLEX, L88-LED-M, PRIMA-FIX, PRIMA-FLEX**

Model code: **ALFA-FIX, ALFA-FLEX, L88-LED-M, PRIMA-FIX, PRIMA-FLEX**

Intended use: SCIALYTIC DIAGNOSTIC LAMP

Risk class of the device, in accordance with the rules laid down in Annex VIII of REGULATION (EU) 2017/745 :

CLASS I

Motivation : Duration: Short term (Annex VIII, CHAPTER I, point 1. DURATION OF USE)
 Description : Non-invasive medical device (Annex VIII, CHAPTER III, point 4 NON INVASIVE DEVICES, sub-section 4.1 Rule 1) Active medical device (Annex VIII, CHAPTER III, point 6. ACTIVE DEVICES, sub-section 6.2 Rule 10)

The manufacturer declares that the device complies with REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 5 April 2017 concerning medical devices and amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC and with the following standards:

- IEC 60601-1 (Part 1: General requirements for basic safety and essential performance)

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 7 on 27	

- IEC 60601-1-2 (Part 2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests)
- IEC 60601-2-41 (Part 1: Particular requirements for the basic safety and essential performance of surgical and diagnostic lighting)

The conformity assessment procedure for the device shall be carried out in accordance with the premises (60) and Article 52 of REGULATION (EU) 2017/745.

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1.5 Warranty certificate

1. The product, including its electrical parts, is guaranteed for 18 months.
2. The warranty takes effect from the date of dispatch of the Product from Verre et Quartz Technologies to the purchaser.
3. In the event of a dispute, the date shown on the 'Transport Document' accompanying the goods will be deemed authentic.
4. The warranty is limited to the supply to the purchaser of replacement parts for the Product or, if Verre et Quartz Technologies considers that the replacement of parts is not feasible, to the replacement of the entire Product; furthermore, it applies to proven manufacturing defects and in the final judgement of Verre et Quartz Technologies. The warranty therefore does not include any other costs or expenses (such as, by way of example and not exhaustively, labour costs, packaging and transport costs, etc.).
5. The guarantee excludes parts subject to natural wear and tear (for example, but not exhaustively: halogen bulbs, LEDs, fuses, relays, ball bearings, etc.).
6. The guarantee does not cover:
 - Malfunctions caused by failure to comply with all the instructions in the manuals;
 - Malfunctions due to installation and/or maintenance errors;
 - Malfunctions or defects due to negligence, improper use or other causes not attributable to Verre et Quartz Technologies;

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 8 on 27	

Malfunctions or defects due to the fact that the electrical installation in the room where the installation takes place does not comply with standard IEC 60364-7-710 (standard for electrical installations in rooms used for medical purposes) and similar standards.

1. Verre et Quartz Technologies will reimburse the purchaser for direct damages documented as being attributable to its Product, caused within the warranty period, for an amount not exceeding 40% of the net value of the Product as shown on the purchaser's invoice. Verre et Quartz Technologies expressly disclaims any liability for indirect or consequential damages (even in cases of non-use of the lamp) arising from the supply.
2. The warranty set out in this certificate replaces statutory warranties for defects and non-conformities and excludes any other possible liability of Verre et Quartz Technologies arising from the Products supplied.
3. Reimbursement for any damage caused to persons or things, due to the malfunction or defects of the Product, will be limited to the limit of Verre et Quartz Technologies' insurance cover for civil liability.
4. The purchaser will automatically forfeit the warranty:
 - if the Product is found to have been altered or modified by the purchaser or third parties;
 - if the Product has been repaired by the purchaser or by third parties without complying the information contained in the Instruction Manuals;
 - if the Product registration number has been erased, hidden or removed;
 - if the purchaser is in arrears with payments.

For repairs under warranty, the purchaser must contact Verre et Quartz Technologies exclusively.

1. Parts replaced under warranty must be returned to Verre et Quartz Technologies, and only if requested by Verre et Quartz Technologies, carriage paid and properly packaged.
2. Any part not returned at the request of Verre et Quartz Technologies will be charged the cost of the part.
3. Verre et Quartz Technologies will not accept parts returned by end users or in any case by parties other than the purchaser.
4. Products returned to Verre et Quartz Technologies must be accompanied by return authorization documentation and a description of the fault.
5. For all matters not covered by this warranty certificate, please refer to Italian law.
6. In the event of any dispute arising out of or in connection with the orders to which this warranty certificate applies and which the parties are unable to resolve amicably, the Court of Milan shall have exclusive jurisdiction.

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 9 on 27	

2 The importance of personal safety

2.1 Purpose of use

The Product is intended to illuminate the area of the patient under observation and diagnosis, and is intended for use in medical clinics.

The Product correctly illuminates the working area from a minimum distance of 40 cm to a maximum distance of approximately 70 cm from the light emission point.

The Product is defined as a diagnostic lamp according to IEC 60601-2-41 :

- A diagnostic lamp is defined as a lamp used for local illumination of the patient's body, with the aim of facilitating diagnosis or treatment which can be interrupted without danger to the PATIENT, in the event of a lack of light. (Not for use in operating theatres).

2.2 Safety conditions (side effects)

Do not point the light source at the eyes of the patient or operator.

- When the use of the product is limited to the face, the patient's eyes must be covered with adequate protection.

Failure to observe these precautions may result in glare and potential damage to the retina.

- Do not place or hang objects on the Product.

Failure to observe this precaution will not guarantee safe positioning and may result in objects falling into the area of use.

- Do not hold on to the Product with all your weight.

Failure to do so may damage the mechanics of the Product.

- Do not cover the dome of the Product during operation.

Failure to do so may prevent heat exchange with the environment and cause the product to overheat.

- Avoid shocks to the moving arms and the head of the Product.

A violent shock could damage the Product and pieces of paint could flake off and fall into the patient's area.

- To avoid any significant risk of reciprocal interference due to the presence of the Product during specific examinations or treatments, refer to section 9.

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 10 on 27	

Power supply unit

To reduce the risk of burns, fire, electric shock or injury to persons or animals:

- Use this power supply only for its intended purpose, as explained later.
- Do not use the mains adapter outdoors; it is designed for indoor use only.
- Do not allow the power pack to be used as a toy. Be careful when the mains adapter is used by or near children.
- Only use accessories recommended by the manufacturer.
- Never use this power supply if its plug or cable is damaged, if it has been dropped or damaged, or if it has been dropped into water. In these circumstances, return the power supply to an authorized dealer or contact customer service.
- Never insert any objects into the openings.
- Never use the power pack in a room where aerosols (sprays) are being used or where oxygen is being administered.
- The power pack should be used close to a convenient and easily accessible mains socket.
- Always disconnect the power pack from the mains immediately after use.

2.3 Environmental conditions

- The product is not intended for use in areas at risk of explosion.
- The product is not intended for use in the presence of flammable mixtures of anaesthetics with air, oxygen or N₂O (laughing gas).
- The Product is not intended for use in an oxygen-rich environment or in the presence of flammable agents.
- During operation, the ambient temperature must be between 10°C and 40°C.
- Relative humidity must be between 30% and 75%.
- Atmospheric pressure must be between 700 and 1060hPa.

2.4 Checks before each use

Before each use, in order to guarantee the safety of the Product and correct diagnosis, the operator must:

- Clean/disinfect the Product in accordance with the provisions laid down by the competent national commission;
- Check that the light emitted is stable and of adequate intensity;
- Check that the mobile arm maintains its position correctly, without falling.

	<h1>USER MANUAL</h1>	12/12/24	<h1>FR</h1>
		Page 11 on 27	

3 Installing the product

Before proceeding with the installation phase, check that all the packaging is present and in good condition, with no transport damage, and that the contents correspond to those listed above.

Complaints will only be considered if the seller or the carrier is informed immediately. All complaints must be made in writing. Goods are always transported at the buyer's risk. Keep the original packaging in case you need to return the Product.

The Product is supplied with a choice of different support systems:

- Table vice (13mm) ;
- Wall bracket (13mm);
- Universal rail clamp (diameter 13mm);
- Heavy rolling stand with 800mm tube (for LED 9)

For the LED 9 model, the pack also includes a sterilisable handpiece.

Do not position the appliance in such a way that it is difficult to reach and disconnect in an emergency.

To avoid any risk of electric shock, this appliance must only be connected to an earthed power supply.

3

3.1 Table-top installation

- Secure the bracket to the table by tightening the threaded pin.
- Insert the lamp into the hole in the top of the bracket.
- Using a screwdriver, tighten the screw at the back of the clamp.

3.2 Wall installation

- Fix the bracket to the wall with 3 expansion screws. Verre et Quartz Technologies does not supply screws.
- The wall must be load-bearing and made of solid bricks. Installation on perforated brick and plasterboard walls is only permitted with the application of a plate on the opposite side of the wall (sandwich closure). Verre et Quartz Technologies suggests using M5 screws.
- Insert the lamp into the hole in the top of the bracket.
- Tighten the threaded knob, ensuring that it fits into the countersink in the lamp stem to prevent accidental release.

3.3 Rail clamp installation

- Attach the rail bar.
- Insert the clamp onto the bar and tighten the lower knob.
- Insert the lamp into the hole in the clamp.
- Tighten the threaded knob, ensuring that it fits into the countersink in the lamp stem to prevent accidental tripping.

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